



# ALL-SILICONE FOLEY CATHETER

- **All-Silicone Construction**  
Soft, comfortable, toxin-free silicone elastomer.
- **Latex-Free**  
Eliminates potential risks that may be associated with repeated exposure to latex devices.
- **Smooth Integrated Balloon**  
The transition between shaft and balloon is smooth and progressive for minimal trauma on insertion and removal.
- **Clear Silicone Shaft**  
Enables inspection of the lumen.
- **Large Diameter Drainage Lumen**  
Improved drainage compared to latex Foley catheters.



# ALL-SILICONE FOLEY CATHETER

The All-Silicone Foley Catheter is designed to enhance comfort and eliminate risks that may be associated with repeated exposure to latex devices. The clear silicone material enables visualization of drainage and helps minimize encrustation. The catheter provides improved drainage with a lumen area that is approximately 30% larger than latex catheters.

## All-Silicone 2-Way Foley Catheter - Paediatric

Individually packaged in 12-pack carton.

BALLOON SIZE	DIAMETER	REF#	PIP CODE
1.5ml	6 Ch	14206	
3ml	8 Ch	14208	
3ml	10 Ch	14210	

## All-Silicone 2-Way Foley Catheter - Standard

Individually packaged in 12-pack carton.

BALLOON SIZE	DIAMETER	REF#	PIP CODE
10ml	12 Ch	14212	301-4628
10ml	14 Ch	14214	301-4636
10ml	16 Ch	14216	301-4644
10ml	18 Ch	14218	301-4651
10ml	20 Ch	14220	301-4669
10ml	22 Ch	14222	301-4677
10ml	24 Ch	14224	301-4685
30ml	16 Ch	24216	301-4560
30ml	18 Ch	24218	301-4578
30ml	20 Ch	24220	301-4586
30ml	22 Ch	24222	301-4594
30ml	24 Ch	24224	301-4602
30ml	26 Ch	24226	301-4610



### DISTRIBUTED BY

Rochester Medical Ltd  
10 Commerce Way  
Lancing, West Sussex BN15 8TA

Freephone: 0800 0322755  
Phone: +44 (0)1903 875055  
Fax: +44 (0)1903 875085

### MANUFACTURED BY

Rochester Medical Corporation, Stewartville, MN 55976 USA

Rochester Medical is a registered trademark of Rochester Medical Corporation.  
US Patent No. 5,137,671.