



MORE CHOICES. MORE ADVANTAGES.



Clear Advantage™ with Aloe




Introducing Clear Advantage™ with Aloe

Self-Adhesive Silicone Sheaths

New Clear Advantage with Aloe offers all the trusted benefits of the original Clear Advantage sheaths, plus a softer, more flexible silicone in three sheath styles. This new generation of products brings security and peace of mind within reach for many men, including those who were previously hard-to-fit.

Silicone Catheters Proven to Breathe Better

Rochester Medical's patented latex-free, clear silicone sheaths allow gentle transpiration and permit skin monitoring while the sheath is in place. By releasing perspiration through the sheath material, silicone sheaths reduce moisture buildup that could lead to maceration and make the sheath more comfortable to wear.



Testing conducted by an independent test lab for Rochester Medical Corporation.

Innovations for a More Comfortable and Secure Fit

Softer Silicone

Like the original design, Clear Advantage with Aloe sheaths continue to reduce the risk of allergic reactions and help ensure optimal skin integrity. The softer silicone also conforms more closely to the contours of the penis, providing better contact for the adhesive and better security against accidental detachment.

More Styles, Each in Five Sizes

There are now three styles of Clear Advantage with Aloe to widen the choice and fit more men – Standard, Pop-On®, and WideBand®.


- The new shorter length Pop-On sheath with forward placement of the adhesive makes it possible to fit men who have a short or retracted penis.
- The new WideBand sheath offers a standard length and a wider band of adhesive, providing reassurance for patients who have concerns about security of attachment.



Style 1 Standard



Style 2 Pop-On



Style 3 WideBand

For more information or to contact your sales representative call Freephone 0800 0322755.



DISTRIBUTED BY Rochester Medical Ltd
10 Commerce Way, Lancing, West Sussex BN15 8TA
Phone: +44 (0)1903 875055 Fax: +44 (0)1903 875085

MANUFACTURED BY Rochester Medical Corporation
Stewartville, MN 55976 USA

Clear Advantage is a trademark of Rochester Medical Corporation. Pop-On, WideBand and Rochester Medical are registered trademarks of Rochester Medical Corporation. Patents: US 5,176,666; US 5,334,175; US 5,376,085; AU 659117; DE 69107470.4; EP 0556325; FR 0556325; GB 0556325; JP 3,409,070.

Introducing Clear Advantage™ with Aloe


Self-Adhesive Silicone Sheaths



New Clear Advantage with Aloe offers all the trusted benefits of the original Clear Advantage sheaths, plus a softer, more flexible silicone in three sheath styles. This new generation of products brings security and peace of mind within reach for many men, including those who were previously hard-to-fit.

Silicone Catheters Proven to Breathe Better

Rochester Medical's patented latex-free, clear silicone sheaths allow gentle transpiration and permit skin monitoring while the sheath is in place. By releasing perspiration through the sheath material, silicone sheaths reduce moisture buildup that could lead to maceration and make the sheath more comfortable to wear.



Testing conducted by an independent test lab for Rochester Medical Corporation.

Innovations for a More Comfortable and Secure Fit

Softer Silicone

Like the original design, Clear Advantage with Aloe sheaths continue to reduce the risk of allergic reactions and help ensure optimal skin integrity. The softer silicone also conforms more closely to the contours of the penis, providing better contact for the adhesive and better security against accidental detachment.

More Styles, Each in Five Sizes

There are now three styles of Clear Advantage with Aloe to widen the choice and fit more men – Standard, Pop-On®, and WideBand®.


- The new shorter length Pop-On sheath with forward placement of the adhesive makes it possible to fit men who have a short or retracted penis.
- The new WideBand sheath offers a standard length and a wider band of adhesive, providing reassurance for patients who have concerns about security of attachment.



Style 1 Standard



Style 2 Pop-On




Style 3 WideBand

For more information or to contact your sales representative call Freephone 0800 0322755.


Clear Advantage™ (original) - Standard

Self-adhesive, clear silicone sheath that is breathable and odourless. Our most popular sheath.

	Size		FP10 PRESCRIPTION CODES	
			30-pack REF #	PIP CODE
	Small	23mm	786187	224-5629
Medium	28mm	786195	224-5637	
Intermediate	31mm	786225	224-5645	
Large	35mm	786233	224-5652	
X-Large	40mm	786241	224-5660	


Clear Advantage™ with Aloe - Style 1, Standard

Self-adhesive, clear silicone sheath with aloe that is soft, breathable and odourless.

	Size		FP10 PRESCRIPTION CODES	
			30-pack REF #	PIP CODE
	Small	24mm	1243	326-2722
Medium	28mm	1283	326-2730	
Intermediate	32mm	1323	326-2748	
Large	36mm	1363	326-2755	
X-Large	40mm	1403	326-2763	


Clear Advantage™ with Aloe - Style 2, Shorter Pop-On®

Self-Adhesive, clear silicone sheath with aloe offers maximum adhesive on a shorter 38mm (1 ½ inch) sheath. Soft, breathable and odourless, the sheath is designed to fit a wider range of patients including those with a short or retracted penis.

	Size		FP10 PRESCRIPTION CODES	
			30-pack REF #	PIP CODE
	Small	24mm	2243	326-2771
Medium	28mm	2283	326-2789	
Intermediate	32mm	2323	326-2797	
Large	36mm	2363	326-2805	
X-Large	40mm	2403	326-2813	

Clear Advantage™ with Aloe - Style 3, Standard With WideBand®

Self-adhesive, clear silicone sheath with aloe offers a 76mm (3 inch) width of adhesive built into a soft, breathable and odourless sheath designed to stay securely in place.

	Size		FP10 PRESCRIPTION CODES	
			30-pack REF #	PIP CODE
	Small	24mm	3243	326-2821
Medium	28mm	3283	326-2839	
Intermediate	32mm	3323	326-2847	
Large	36mm	3363	326-2854	
X-Large	40mm	3403	326-2862	



DISTRIBUTED BY

Rochester Medical Ltd
10 Commerce Way
Lancing, West Sussex BN15 8TA

Freephone: 0800 0322755
Phone: +44 (0)1903 875055
Fax: +44 (0)1903 875085

MANUFACTURED BY

Rochester Medical Corporation, Stewartville, MN 55976 USA